MANUEL D'UTILISATION • USER GUIDE • GEBRAUCHSANLEITUNG • MANUAL DE USO • ISTRUZIONI D'USO • INSTRUÇÕES DE USO  
BRUGERVEJLEDNING • BRUKSANVISNING • GEBRUIKSAANWIJZING • INSTRUKCJA OBSŁUGI • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ  
ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ • UŽIVATELSKÝ MANUÁL • 使用说明 • 使用説明書


# Respireo SOFT Baby

Masque nasal avec fuites calibrées • Nasal mask vented • Nasenmaske mit kalibrierten luftauslassöffnungen • Mascarilla nasal con orificios de exhalación calibrados • Maschera nasale con fori di esalazione calibrati • Máscara nasal com furos de exalação calibrados  
Nasal maske med kalibrerede udåndingshuller • Näsmask med kalibrerade utandningshål • Nesemaske med kalibrerte utpustingshull  
• Neusmasker met gekalibreerde ventilatieopeningen • Maska nosowa z kalibrowanymi otworami wentylacyjnymi • Ρινική μάσκα με ρυθμιζόμενες οπές εξαερισμού • Носовая маска с калиброванными выдыхательными отверстиями • Nosní maska s kalibrovanými vydechovacími otvory • 带有标准通气口的鼻面罩 • 鼻マスク 通気穴付き


FRANÇAIS.....	p. 7
ENGLISH.....	p. 12
DEUTSCH.....	p. 16
ESPAÑOL.....	p. 21
ITALIANO.....	p. 26
PORTUGUÊS.....	p. 31
DANSK.....	p.37
SVENSKA.....	p. 41
NORSK.....	p. 45
NEDERLANDS.....	p. 50
POLSKI.....	p. 55
ΕΛΛΗΝΙΚΑ.....	p. 61
РУССКИЙ.....	p. 66
ČEŠTINA.....	p. 70
中文.....	p. 74
日本語.....	p. 79

# Respireo SOFT Baby


**OPTIONAL**


1


2


3


4



2


3


4

1

2

**1 - USAGE PREVU**

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## 1 - APPLICATION

Respireo SOFT Baby is a non-invasive nasal mask used for treating respiratory failure in babies weighing between 3 and 12 kg.

Respireo SOFT Baby is designed for use by babies who have been prescribed continuous positive airway pressure (CPAP) or Bi-level positive airway pressure treatment.

Respireo SOFT Baby is designed for use by a single patient at home or more than one in the hospital environment.

## ⚠ 2 - IMPORTANT NOTES

- 2.1 The mask must be used under the supervision of an adult who has read and understands the instructions.
- 2.2 The mask must only be used with a CPAP device or Bi-level flow generator recommended by the family doctor or respiratory therapist, operating at a pressure of at least 4 cmH<sub>2</sub>O. With low CPAP pressures, the flow through the vent ports may be insufficient to eliminate all the gas exhaled from the tube, which could lead to partial re-inhalation.
- 2.3 The vents in the mask must never be obstructed.
- 2.4 The operating pressure range is 4-24 cmH<sub>2</sub>O.
- 2.5 If you experience discomfort or irritation, or suffer an allergic reaction to any of the components of the mask, please report this to your doctor or respiratory therapist.
- 2.6 If oxygen is used with the device (CPAP or Bi-level), the oxygen supply must be turned off when it is not in use to prevent oxygen accumulating inside, which is a fire risk.
- 2.7 Never smoke or light naked flames when using oxygen.
- 2.8 Do not use the mask in the presence of nausea or vomiting.

- 2.9 The mask must be cleaned and disinfected. For further details refer to section CLEANING THE MASK.
- 2.10 Check the mask each time before use. If there is any visible damage, report it to the dealer.
- 2.11 The mask must only be used in accordance with these instructions only. No other instructions are applicable.
- 2.12 Some of the components of the mask are small enough to be swallowed by children and must never be left unattended. An adult must always be present.

## 3 - COMPONENTS OF THE MASK (Fig.A)

- (a) Main body:** frame (1), cushion (2), forehead rest (3);  
**(b) Tube assembly:** elbow fitting with calibrated vent ports (4), tube (5), swivel fitting (6), cap (7);  
**(c) Headgear.**

## 4 - HOW TO ASSEMBLE THE MASK (Fig. B-C)

Connect the tube assembly (b) by inserting the sphere on the elbow fitting (4) in the slot in the main body (a), as shown in Fig. B1.

The instructions on how to use the mask depend on the chosen configuration (WITH or WITHOUT the forehead rest). Refer to the relevant subsection:

- 4.1 – Using the mask WITH the forehead rest (3)
- 4.2 – Using the mask WITHOUT the forehead rest (3)

### 4.1 USING THE MASK WITH THE FOREHEAD REST (FIG. B2)

4.1.1 Insert the three straps at the top of the headgear (c) in the three slots in the forehead rest (3), as shown in Fig. B2. Then fasten the straps.

4.1.2 Move on to Section HOW TO PUT THE MASK ON THE CHILD.

## **4.2 USING THE MASK WITHOUT THE FOREHEAD REST (FIG. C) (OPTION)**

4.2.1 To remove the forehead rest (3), merely cut it off the main body (a) using a pair of scissors, as shown in Fig. C1.

*N.B. Do not throw away the forehead rest (3). It will be required later.*

*N.B. Once the forehead rest (3) has been cut off from the main body (a), it cannot be reattached.*

4.2.2 Take the cut-off forehead rest (3) and insert the central strap at the top of the headgear (c) in the top slot, as shown in Fig. C2.

4.2.3 Insert the side straps at the top of the headgear (c) in the side slots, as shown in Fig. C3.

4.2.4 Insert the central strap at the top of the headgear (c) in the slot at the top of the main body (a), as shown in Fig. C4. Then fasten the strap.

## **5 - HOW TO PUT THE MASK ON THE CHILD (Fig.D)**

5.1 Place the mask gently on the child's face, making sure the cushion (2) fits snugly round the nose. Pull the headgear over the child's head, as shown in Fig. D1.

5.2 Insert the loose lower straps of the headgear (c) in the side slots in the main body (a), as shown in Fig. D2. Then tighten the straps.

5.3 Position the headgear (c) comfortably on the child's head and adjust all the straps, starting from the lower ones, as shown in Fig. D3.

5.4 Connect the swivel fitting (6) to the tube of the CPAP/Bi-level device, as shown in Fig. D4. Then turn on the flow generator.

5.5 Readjust the lower straps if necessary to eliminate any leaks, but make sure the mask is not too tight.

*N.B. The lower straps only serve to hold the mask in position, so they must be kept fairly loose.*

## **6 - HOW TO REMOVE THE MASK**

Unfasten one of the lower straps on the headgear (c). Take off the mask off the face and lift the headgear off. It is not necessary to unfasten the other straps.

## **7 - HOW TO DISSAMBLE THE MASK (Fig. E)**

7.1 Unfasten and pull out the straps, as shown in Fig. E1, then detach the headgear (c) from the main body (a).

7.2 Detach the tube assembly (b) by pulling the elbow fitting (4) out of the main body (a), as shown in Fig. E2.

## **8 - HOW TO CLEAN THE MASK**

### **8.1 – DAILY CLEANING FOR REUSE BY THE SAME PATIENT**

First wash your hands thoroughly.

Disassemble all the components (see HOW TO DISASSEMBLE THE MASK) and wash carefully in warm water at about 30°C, using mild soap. Rinse thoroughly under the tap and leave in the air to dry, but not in direct sunlight.

*N.B. The headgear (c) should be cleaned once a week.*

### **8.2 – REPROCESSING THE MASK FOR RE-USE WITH OTHER PATIENTS**

*N.B. The mask must always be regenerated before being used by another patient.*

*N.B. The headgear (c) must be replaced each time.*

*N.B. In case of thermal disinfection, chemical disinfection and autoclave sterilization, the tube assembly (b) must be replaced.*

The mask can be regenerated up to 12 times, using one of the following methods:

- high-level thermal disinfection;
- high-level chemical disinfection;
- sterilisation in an autoclave at 121°C for no more than 20 minutes;
- Sterrad gas-plasma sterilisation.

Full details of the mask regeneration process can be obtained by logging on to [www.device.airliquidehealthcare.com](http://www.device.airliquidehealthcare.com) or by contacting your dealer.

### 9 - PRECAUTIONS

Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents or aromatic oils as they would damage the product and reduce its life.

Prolonged exposure to direct sunlight or heat may damage the mask.

If any damage is noticed (cracks, tears, etc.), the damaged part must be replaced.

Check carefully to make sure there are no soap deposits or foreign bodies inside the mask.

### 10 - PRESSURE/OXYGEN PORT

If additional oxygen or a pressure reading is required, use the port in the elbow fitting (4).

Remove the cap (7) and insert the oxygen/pressure tube in the port.

## 11 - TROUBLE-SHOOTING CHART

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
<ul style="list-style-type: none"> <li>• No air or insufficient air comes out of the mask</li> </ul>	<ul style="list-style-type: none"> <li>• Flow generator off</li> <li>• Flow generator regulated wrongly</li> <li>• Patient tube not connected</li> </ul>	<ul style="list-style-type: none"> <li>• Switch on the flow generator.</li> <li>• Regulate the flow generator.</li> <li>• Connect the patient tube to the flow generator and mask.</li> </ul>
<ul style="list-style-type: none"> <li>• The mask is noisy</li> </ul>	<ul style="list-style-type: none"> <li>• Vent ports partially or entirely obstructed</li> <li>• Mask assembled incorrectly</li> </ul>	<ul style="list-style-type: none"> <li>• Clean the vents using a SOFT brush.</li> <li>• Take the mask apart and reassemble it, as explained in section HOW TO ASSEMBLE THE MASK.</li> </ul>
<ul style="list-style-type: none"> <li>• Exhaled air does not flow out and difficulty exhaling</li> </ul>	<ul style="list-style-type: none"> <li>• Vent ports partially or entirely obstructed</li> </ul>	<ul style="list-style-type: none"> <li>• Clean the vents using a SOFT brush.</li> </ul>
<ul style="list-style-type: none"> <li>• The mask does not fit the face snugly</li> </ul>	<ul style="list-style-type: none"> <li>• Mask not positioned correctly on the patient's face</li> <li>• Headgear straps too tight</li> </ul>	<ul style="list-style-type: none"> <li>• Reposition the mask, as explained in section HOW TO PUT THE MASK ON THE CHILD.</li> <li>• Unfasten the lower straps and refasten more loosely.</li> </ul>


## 12 - DISPOSAL INSTRUCTIONS

The components of the mask are made of the following materials and must be disposed of in accordance with the applicable laws.

COMPONENT	MATERIAL
Main body (a)	Silicone + Polycarbonate
Tube assembly (b)	Polycarbonate + Metallocene + Acetalic + Silicone
Headgear (c)	Nylon + polyurethane

## 13 - TECHNICAL DATA

### MINIMUM CONTINUOUS FLOW

*N.B. Due to manufacturing differences, the exhalation flow rate may differ from the values shown in the table ( $\pm 10\%$ ).*

Normative references	
	CE marking in accordance with European Directive 93/42/EEC, class IIa, as amended, EN ISO 14971, ISO 5356-1, ISO 10993-1 and EN ISO 17510-2, EN ISO 17664

<b>Resistance</b>	1,05 cmH <sub>2</sub> O at 50 lpm – 3,30 cmH <sub>2</sub> O at 100 lpm
<b>Dead space volume inside the mask</b>	63 ml
<b>Weighted sound pressure A according to ISO 3744</b>	15,8 dB(A) $\pm$ 1,5 dB(A)
<b>Weighted acoustic pressure A according to ISO 3744</b>	19,5 dB(A) $\pm$ 1,5 dB(A)
<b>Operating pressure</b>	4-24 cmH <sub>2</sub> O
<b>Temperature range</b>	5-40°C
<b>Max. overall dimensions</b>	86 mm x 70 mm x 352 mm

## 14 - WARRANTY

Air Liquide Medical Systems S.r.l. warrants the product free from material and manufacturing defects for a period of 120 days from the date of purchase, provided the operating conditions and instructions are complied with.

If the product is found to be faulty during normal use, it will be repaired or replaced at the manufacturer's discretion.

The manufacturer accepts liability for the safety, reliability and performance of the device only if it is used for the intended application and in accordance with the operating instructions.

The manufacturer accepts liability for the safety, reliability and performance of the device only if repairs are carried out by the manufacturer or an authorised workshop.

The original packaging must be retained for the entire warranty period.

Claims under the warranty should be submitted to your local dealer, who can also provide further details of your warranty rights.

**COPYRIGHT**

The information contained in this booklet must not be used for any purpose other than the reason for which it is provided. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

**TECHNICAL UPDATES**

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated to make sure they remain consistent with the features of the devices placed on the market. If the booklet accompanying the device gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

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AIR LIQUIDE MEDICAL SYSTEMS S.p.A.