
DE Gebrauchsanweisung EN Instructions for Use FR Mode d'emploi
NL Gebruiksaanwijzing IT Istruzioni d'uso TR Kullanma Kılavuzu
ES Manual de instrucciones

WM 67480F 09/2016 DE, EN, FR, NL, IT, TR, ES(MX)







CE 0197

Löwenstein Medical
Technology GmbH + Co. KG
Kronsaalweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.de

JOYCE SilkGel

Nasal Mask
vented, non vented

1 Bedienung

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7 Störungen

Table with 3 columns: Störung, Ursache, Behebung. Lists common issues like mask fit, air leakage, and pressure measurement.

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Table with 2 columns: Parameter, Value. Lists technical specifications like dimensions, weight, and pressure range.


Table with 2 columns: Parameter, Value. Lists flow resistance and noise emission data.

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1 Operation

The following figures show you how to adjust, apply, remove, dismantle and assemble the mask:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask
- 4 Dismantling the mask
- 5 Assembling the mask

2 Introduction

2.1 Intended use

The JOYCE SilkGel mask is used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device.

2.2 Contraindications

The mask may not be used, or may be used only with particular caution, if any of the following symptoms are present:

erosions and ulcerations, skin allergies, flush of the facial skin, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, intake of medication which may induce vomiting, the necessity for immediate intubation. Please also observe the contraindications in the instructions for use of your therapy device.

2.3 Side effects

Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

3 Safety

Risk of injury from the supply of oxygen!

Oxygen can become deposited in clothing, bedlinen and hair. In conjunction with smoking, naked flames or electrical equipment, it can cause fires and explosions.

- ⇒ Do not smoke.
- ⇒ Avoid naked flames.

Risk of injury from re-inhaling CO₂!

If the mask is used incorrectly, CO₂ may be re-inhaled.

- ⇒ Do not close off the exhalation system of the mask.

- ⇒ Only put on the mask for an extended period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who are unable to take it off themselves.

Risk of injury if the mask slips!

- If the mask slips or falls off, the therapy is ineffective.
- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.

Risk of injury from anesthetic gases!

- Anesthetic gas may escape through the exhalation valve and put third parties at risk.
- ⇒ Never use the mask during anesthesia.

4 Product description

4.1 Overview

A diagram of the individual parts can be found on the title page.

1. Headgear
2. Forehead cushion
3. Forehead support
4. Retaining ring
5. Elbow
6. Pressure measurement port seal
7. Pressure measurement port
8. Rotating sleeve
9. Mask body
10. Headgear clip
11. Mask cushion
12. Tube anchoring strap
- 13.

4.2 Compatible devices

With many combinations of device, actual pressure in the mask will not correspond to prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that actual pressure in the mask corresponds to therapy pressure.

4.3 Exhalation system

vented variants

Vented variants (transparent elbow and rotating sleeve) have an integrated exhalation system. The retaining ring and mask body are shaped so that there is a gap between these parts. The exhaled air can escape through this gap.

Non-vented variants

Non-vented variants (blue retaining ring and elbow) do not have an exhalation system. If you are using the non-vented variant with an external exhalation valve, please see the associated instructions for use.

4.4 Pressure measurement port

The mask has a pressure measurement port, which can be used for measuring pressure or supplying oxygen. If you are not using the connection, close it with the seal for pressure measurement to reach therapy pressure.

5 Hygiene treatment

⚠ WARNING

Risk of injury from inadequate cleaning!

Residues may congest the mask, impair the integrated exhalation system and jeopardize therapy success.

⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

5.1 Clean mask

1. Dismantle mask (see Figure 4).
2. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	X	
When washing mask parts (except for mask cushion and forehead cushion), clean thoroughly with a cloth or a soft brush.		X
Wash headgear by hand.		X

3. Rinse all parts with clear water.
4. Allow all parts to air-dry.
5. Perform a visual inspection.
6. If necessary: replace damaged parts.

7. Re-assemble mask (see Figure 5).

i	Discolorations of mask components do not impair the functionality of the mask.
----------	--

5.2 Change of patient

In the clinical sphere only: You must subject the mask to a hygiene treatment in the event of a change of patient. Information on the hygiene treatment with a change of patient can be found in a brochure available on the manufacturer's website. We will also send you this brochure on request.

6 Disposal

You can dispose of all parts in domestic waste.

7 Troubleshooting

Fault	Cause	Remedy
Pain from pressure on the face.	Mask too tight.	Loosen headgear slightly.
Draft in the eye.	Mask too loose.	Tighten headgear slightly.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure is not reached.	Mask not correctly adjusted.	Re-adjust mask (see Figure 2).
	Mask cushion amaged.	Replace mask cushion.
	Patient circuit leaking.	Check connectors, check tubes properly located.
	Air escapes at pressure measurement port.	Close pressure measurement port with plug (included in scope of supply).

8 Technical specifications

Product class to directive 93/42/EEC	Ila
Dimensions (W x H x D)	80 mm x 130 mm x 90 mm
Weight	85 g
Therapy pressure range	4 hPa - 40 hPa
Tube connection: tapered connection to EN ISO 5356-1 vented non vented	Ø 22 mm (male) Ø 22 mm (female)
Pressure measurement port	Ø 4 mm


Temperature range: Operation Storage	+5 °C to + 40 °C -20 °C to +70 °C
Flow resistance vented at 50 l/min at 100 l/min	0,02 hPa 0,26 hPa
Flow resistance non vented at 50 l/min at 100 l/min	0,09 hPa 0,37 hPa
Quoted two-figure noise emission value to ISO 4871: - Sound pressure level - Sound power level - Uncertainty factor	15 dB(A) 23 dB(A) 3 dB(A)
Service life	up to 12 months ¹
Standards applied	EN ISO 17510-2:2009

¹ The mask materials deteriorate if exposed to e.g. aggressive detergents. In individual cases it may be necessary to replace mask parts sooner.

The right to make design modifications is reserved.

9 Pressure/flow curve

The characteristic pressure/flow curve shows the leakage flow as a function of therapy pressure.

10 Materials

Mask part	Material
Mask cushion, forehead cushion	Silicone, silicone gel, silicone coating
Plug for pressure measurement port	Silicone
Retaining ring	TPE (thermoplastic elastomer), PP (polypropylene)
Mask body, elbow	PC (polycarbonate), PA (polyamide)

Mask part	Material
Headgear clip	POM (polyoxymethylene)
Forehead support, rotating sleeve	PA (polyamide)
Headgear	Elasthane, polyester, PU (polyurethane), UBL loop, cotton, PA (polyamide)

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

11 Warranty

Löwenstein Medical gives the customer a limited manufacturer warranty on new original Löwenstein Medical products and any replacement part fitted by Löwenstein Medical in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions are available on the website of the manufacturer. We will also send you the warranty conditions on request.

Please bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine replacement parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

12 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

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