

NEW!

OXYvent Cube 30 ATV

- Non-invasive Ventilator -

**Special price
 on request**

The OXYvent Cube 30 ATV is a device for non-invasive, non-life supporting respiration of spontaneously breathing patients who have a body weight of at least 13 kg. Especially worthy of mention is the unique design and the ease of use of the device - **made in Germany.**

● Indications

- COVID-19
- COPD
- Neurological, muscular and neuromuscular complaints (e.g. paresis of the diaphragm)
- Restrictive ventilation disorders (e.g. scoliosis, thorax deformities)
- Central respiratory regulatory disorders
- Obstructive sleep apnea (OSA)

● Benefits at a glance


- Functional and easy-to-use design
- Optimal overview and high user-friendliness due to the large 4,3" colour display
- High sensitive flow based trigger detection with trigger failure auto-detect function
- Modern target volume function - can be set based on breath-by-breath target volume or target minute volume
- Trigger lock adjustable
- Extensive device alarm management
- Two independent adjustable therapy-settings possible (e.g. therapy mode for day and night)
- Production and development according to the highest German quality standards

● Technical Data

Ventilation modes	CPAP, BiLevel S/PSV-S, BiLevel T/PCV, BiLevel ST/PSV, APCV
Pressure range	3 - 30 mbar
CPAP-pressure range	3 - 20 mbar
Target volume	200 - 2.000 ml
Target minute volume	2.000 - 15.000 ml
Trigger	Insp.: 5 level settings, Exsp.: auto, manual 70-30%
Trigger Lock	0,3 - 8 sec.
Volume protection	Target volume: 200 - 2.000ml, Target minute volume: 2.000 - 15.000ml

● Scope of delivery

- Cube 30 ATV device
- Humidifier (heated)
- 15mm tube system
- EU power cable
- Fine filter
- Coarse dust filter
- Transport bag
- Manual



**Target volume and
 Target minute volume
 adjustable**